



# HORIBA

Explore the future





# Fluorescence in Biomolecular Research

HORIBA Jobin Yvon IBH Ltd  
A division of HORIBA Scientific

Kulwinder Sagoo  
Product Specialist  
February 27, 2013

# Outline



- History of IBH
- HORIBA Group
- What is Fluorescence?
- Applications relevant to Biomolecular research
- Equipment

# Company Timeline

- Founders – David Birch, Bob Imhof and Tony Hallam
- Spin off company from Strathclyde University
- Role of IBH - to design fluorescence lifetime spectrometers

**1977-** IBH Consultants Ltd. formed and designs licensed



**1989-** IBH independent marketing & manufacturing begins




**2003-** IBH joins Jobin Yvon Inc.



**2004-** Name changed to HORIBA Jobin Yvon IBH Ltd.



**2008 -** HORIBA Scientific Group formed



# HORIBA

**HORIBA**  
Scientific



- \$1.5 Billion sales
- Over 5000 employees
- Key business segments:
  - Automotive
  - Semiconductor systems
  - Scientific instruments
  - Environmental analyzers
  - Medical




Head Office – Kyoto, Japan

# HORIBA Scientific

**HORIBA**  
Scientific



## 3 Centers of Excellence:



### **HOR (Kyoto, head office)**

X-ray products, elemental analyzers (C,N,H,O, S), particle analyzers


### **HIC (JYUS):**

Fluorescence (includes IBH in Glasgow), spectroscopic components, forensics, OEM


### **JY SAS:**

Raman, elemental analyzers (ICP, GD), ellipsometry, OEM, SPRi

# Fluorescence Division



2 Centers of Excellence



# What is Fluorescence?

Fluorescence is a multiparameter signal

$$FL = f(I, \lambda_{exc}, \lambda_{em}, p, x, t)$$

$I$  = intensity - measurement is quantum yield ( $\phi$ ),

$\lambda_{exc}$  = excitation wavelength

– measurement of absorption spectrum,

$\lambda_{em}$  = emission wavelength


– measurement of fluorescence spectrum,

$p$  = polarisation – measurement of anisotropy,

$x$  = position – measurement by fluorescence microscopy,


$t$  = time – measurement of fluorescence lifetime.

# Fluorescence




Radiative deactivation of the first excited singlet state

depends on molecule and environment



$$I = I_0 10^{-\epsilon[c]d}$$

$\epsilon[c]d$  = absorbance  
 $\epsilon$  extinction coefficient  
 $[c]$  concentration  
 $d$  pathlength



quantum yield  $\phi$

Stokes shift

Decay time  $\tau$

# Fluorescence techniques

Fluorescence measurement aims to record one or more of these parameters

intensity, wavelength, time, polarisation, position (x,y,z)

## Fluorescence

steady state  
(average)



$$\phi = \frac{\text{emitted photons}}{\text{absorbed photons}} \\ = k_r / (k_r + k_{nr})$$

time-resolved

$$\tau\text{- fluorescence lifetime} \\ = 1 / (k_r + k_{nr})$$



**Absolute measurement**

# The fluorescence lifetime



$$\tau = \frac{1}{k_r + k_{nr}}$$

# Practical set up



[www.picocomponents.com](http://www.picocomponents.com)

# Applications

## ■ Examples

- Fast acquisition of short lifetimes
- FRET (Foerster Resonance Energy Transfer)
- Stern-Volmer quenching
- Lanthanide luminescence
- Time-resolved anisotropy
- Protein Fluorescence
- Solar cell analysis
- Singlet oxygen measurements
- Materials research
- Photophysical research

# Why measure lifetimes?

- Reveal information about kinetic rates, dynamics and structure of fluorophore
- Single molecule sensitivity
- Independent of sample intensity
  - Absolute measurement
  - Concentration independent (within limits)
- Sensitive to local environment
  - pH, local charged group, quenchers
  - Temperature, polarity, viscosity
- Extra Specificity
  - Discriminate against unwanted fluorescence and excitation
- Other information
  - Additional dimension to fluorescence data map
  - Increases specificity of the measurement

# Why measure lifetime in BioScience?

---

- Lifetime time scale of <math><1\text{ns}</math> to   - allows probing of the environment surrounding the fluorophore
    - nm distances (    - ns timescale (
- Also, very sensitive, works in-situ, & non-destructive
- Fluorophores can be intrinsic or extrinsic

# Examples of Fluorescence in Biochemistry & Medicine



- Fluorescence Microscopy
  - Tissues & cells - label antibodies with fluorophore and monitor interaction with antigen
  - FLIM (Fluorescence lifetime imaging microscopy) – bio-molecular interactions resulted in changes in fluorescence lifetime
- FRET – to monitor ie. protein interactions and for biosensors
- Biotechnology – develop biosensors for detection of glucose monitoring
- DNA – sequence sorting by labelling specific terminal bases with a specific fluorophore (chain termination method)

# Questions related to Biomolecular research



How large is the biomolecule of interest?

Has my protein denatured?

Are you exploring new sensitizers for photodynamic therapy?

What proportion of my sample is bound?

Is the biomolecule thermally stable?

Has a protein-protein interaction occurred?

Has the drug bound to the active region of the biomolecule?

Has the molecular microenvironment changed?

How mobile is my biomolecule after a binding interaction has occurred?

Has immobilising my biomolecule affected conformation and function??

Is the active site of the biomolecule accessible for binding ?

Has there been a conformational change?

# How can fluorescence lifetimes help?



## STRUCTURE/ CONFORMATION

### Monitor:

- Changes in conformation because of unfolding, denaturing or binding
- Inter vs Intra molecular distances
- Photosensitisation

### Using techniques such as:

- Fluorescence lifetime
- TRES
- FRET
- TR anisotropy
- Phosphorescence lifetime

## SIZE / MOBILITY

### Monitor:

- Viscosity
- Rotational diffusion
- Binding
- Restricted mobility

### Using techniques such as:

- Fluorescence lifetime
- FRET
- TR Anisotropy

## FUNCTION

### Monitor:


- Binding of proteins, ligands, drugs, etc.
- Bound vs unbound
- Change in rate of rotation of species upon binding

### Using techniques such as:

- Fluorescence lifetime
- TR anisotropy
- Fluorescence quenching

# Measurement Example

## Monitoring viscosity changes



Mitochondrial probe



Stillbenoid molecule (DASPMI)– molecular rotor sensitive to changes in local nanoscale viscosity


Viscosity dependent intermolecular re-arrangement that affects lifetime

G. Hungerford et al. **2009**. *Monitoring sol to gel transitions via fluorescence lifetime determination using viscosity sensitive fluorescent probes*. J. Phys. Chem. B. **113**, 12067-12074.



# Viscosity Study- Cont'd

## ■ Polysaccharides

- Repeating units of monomers linked by glycosidic bonds
- Important role in animal and plant nutrition and structure
- Gellan gum is used in food industry, drug delivery and tissue engineering




DASPMI entrapped in Gellan Gum analysed using FLIM



Viscosity changed ~600x moving from dry to wet front

G. Hungerford et al. **2012**. *In-situ formation of silver nanostructures within a polysaccharide film and its application as a potential biocompatible fluorescence sensing medium*. *Soft Matter*. **8**, 653-659.

# Measurement Example



Picosecond lifetimes


Fast Acquisition

Example of fast data acquisition with DeltaDiode excitation with 100MHz repetition rate and low deadtime (<10ns) timing electronics (Bodipy Derivative)

D. McLoskey et al. 2011. *Fast time-correlated single-photon counting fluorescence lifetime acquisition using a 100 MHz semiconductor excitation source*. Meas. Sci. Technol. 22, 067001.

# Measurement Example- Binding Study

## HSA binding to curcuminoid





### Kinetic TCSPC


- 10 000 decays, 10ms/decay
- $\langle \tau \rangle$  vs macrotime
- $\langle \tau \rangle$  increases upon binding
- total volume - 5 $\mu$ l



- Protein binding important to understand interactions with variety of target molecules
- Phenolic compounds provide anti-oxidant activity with potential health benefits



# Fluorescence Lifetime systems




# Our New Delta Series

**HORIBA**  
Scientific



■ DeltaPro

■ DeltaFlex



[www.deltatcpsc.com](http://www.deltatcpsc.com)

# FluoroFest

## Turn **Fluorescence** into the Sharpest Tool in Your Lab



### Fluorescence Workshop & Course

**Bio and medical research** -  
Nanomedicine, protein, skin, metabolic monitoring, cells, metal ions, glucose sensing

**Nanotechnology** -  
Nanometrology, nanoparticles, quantum dots, carbon nanotubes

**Materials science** -  
Sol-gels, polymers, ceramics, composites, colloids, intelligent and green materials

■ **Learn** the State of the Art in Fluorescence Techniques from world-famous Fluorescence Experts

■ **Experience** hands-on operation of State-of-the-Art Equipment, both lifetime and steady-state, from the supplier of the World's Most Sensitive Spectrofluorometers

■ **Present** your research; exchange ideas with your colleagues.

**Prague, Czech Republic**  
**June 3-5 2009**

For information and registration email:  
Geraldyn.Caruso@jobinyvon.com

*FluoroFest Workshop*

## Turn **Fluorescence** into the workhorse in your lab



### Fluorescence Workshop & Course

**Bio and medical research** -  
Nanomedicine, protein, skin, metabolic monitoring, cells, metal ions, glucose sensing

**Nanotechnology** -  
Nanometrology, nanoparticles, quantum dots, carbon nanotubes

**Materials science** -  
Sol-gels, polymers, ceramics, composites, colloids, intelligent and green materials

■ **Learn** the State of the Art in Fluorescence Techniques from world-famous Fluorescence Experts

■ **Experience** hands-on operation of State-of-the-Art Equipment, both lifetime and steady-state, from the supplier of the World's Most Sensitive Spectrofluorometers

■ **Present** your research; exchange ideas with your colleagues.

**Munich, Germany**  
**Nov 29-Dec 1, 2011**

For information and registration see:  
[fluorofest.org](http://fluorofest.org)

*FluoroFest Workshop*

# Contact

---

**HORIBA**  
Scientific



HORIBA Jobin Yvon IBH Ltd.

Skypark 5, Floor 1

45 Finnieston Street

Glasgow

0141-229-6784

[www.horiba.com](http://www.horiba.com)

